




### Option 1: Share using Dexcom CLARITY app

- 1 Install the Dexcom CLARITY app on your smartphone.
- 2 Log in with your Dexcom account.  
If you use Dexcom Mobile, you'll be logged in automatically.
- 3 Tap Share.
- 4 Tap Accept Invitation.
- 5 Enter the sharing code provided by your clinic.
- 6 Tap Enter Birthday to select your date of birth.
- 7 Tap Continue.
- 8 Tap the I consent to share my data with my clinic box to select it.
- 9 Tap Yes, Share My Data.



### Option 2: Share using CLARITY website

- 1 Receive a printed or emailed sharing invitation from your clinic.
- 2 Enter <https://clarity.dexcom.com/share> in your internet browser and click Log in.
- 3 Enter your Dexcom or Dexcom Mobile login, then click Sign in.
- 4 Enter the invitation sharing code and your date of birth.
- 5 Click I consent to share my data with my clinic.
- 6 Click Yes, Share My Data.



Minimum requirements: [dexcom.com/clarity-requirements](https://dexcom.com/clarity-requirements) | Dexcom Support: 1-888-738-3646, [techsupport@dexcom.com](mailto:techsupport@dexcom.com), or [dexcom.com/support](https://dexcom.com/support)

#### Intended Use/Safety Statement

The web-based Dexcom CLARITY® software is intended for use by both home users and healthcare professionals to assist people with diabetes in the review, analysis, and evaluation of historical CGM data to support effective diabetes management. It is intended for use as an accessory to Dexcom CGM devices with data interface capabilities. Caution: The software does not provide any medical advice and should not be used for that purpose. Home users must consult a healthcare professional before making any medical interpretation and therapy adjustments from the information in the software. Caution: Healthcare professionals should use information in the software in conjunction with other clinical information available to them. Caution: Federal (US) law restricts this device to sale by or on the order of a licensed healthcare professional.