

PRODUCT ALERT

ENFit


OLD

Change From:

**CURRENT SYRINGE,
OR TUBE SET**

AND

**CURRENT
FEEDING TUBE OR
EXTENSION SET**



NEW


Change To:

ENTERAL SPECIFIC

ENFit SYRINGE

AND

**ENFit
FEEDING TUBE OR
EXTENSION SET.**



National Safety organizations are requiring a change in the design for enteral (tube feeding) connectors to improve patient safety. Companies across the industry helped to develop the **ENFit Connection System**. For more information you may visit stayconnected.org

The new **ENFit** connector has a unique design that provides a simple way to reduce the risk of tube feeding misconnections. The **ENFit** design on the ends of tubing, syringes feeding tubes and extension sets allows all the enteral products to fit together, but prevents any other type of connector from fitting with the enteral products.

The nursing staffs in the hospital and homecare will work with you to be sure that you understand this product change, and will be certain to make sure that you have the materials that you need to take care of your child.